

REFERENCE

NEEDLE COMPARISON CHARTS

RAZOR POINT™

PRECISION COSMETIC CONVENTIONAL CUTTING NEEDLES ▲


| NEEDLE | CP MEDICAL | ETHICON | COVIDIEN | SHARPOINT | LOOK |
|-------------------|------------|---------|----------|-----------|------|
| 13mm • 3/8 Circle | NPC-1 | PC-1 | DX-13 | DSCM13 | • |
| 16mm • 3/8 Circle | NPC-3 | PC-3 | DX-16 | DSCM16 | • |
| 19mm • 3/8 Circle | NPC-5 | PC-5 | DX-19 | DSCM19 | • |



RAZOR POINT™


PRECISION REVERSE CUTTING NEEDLES ▼

| NEEDLE | CP MEDICAL | ETHICON | COVIDIEN | SHARPOINT | LOOK |
|-------------------|------------|---------|----------|-----------|------|
| 8mm • 3/8 Circle | NP-6 | P-6 | P-16 | • | • |
| 11mm • 3/8 Circle | NP-1 | P-1 | P-10 | DSM11 | C1 |
| 13mm • 3/8 Circle | NP-3 | P-3 | P-13 | DSM13 | C3 |
| 16mm • 3/8 Circle | NPS-3 | PS-3 | P-11 | DSM16 | • |
| 19mm • 3/8 Circle | NPS-2 | PS-2 | P-12 | DSM19 | PC31 |
| 24mm • 3/8 Circle | NPS-1 | PS-1 | P-14 | DSM24 | PC33 |
| 26mm • 3/8 Circle | NPS | PS | • | • | • |
| 8mm • 1/2 Circle | NP-2 | P-2 | • | • | • |
| 16mm • 1/2 Circle | NPS-4 | PS-4 | P-24 | • | PC30 |



REVERSE CUTTING NEEDLES ▼

| NEEDLE | CP MEDICAL | ETHICON | COVIDIEN | SHARPOINT | LOOK |
|-------------------|------------|---------|----------|-----------|------|
| 13mm • 3/8 Circle | NC-3 | C-3 | C-1 | DS12 | C17 |
| 16mm • 3/8 Circle | NFS-3 | FS-3 | C-12 | DS16 | C22 |
| 19mm • 3/8 Circle | NFS-2 | FS-2 | C-13 | DS18 | C6 |
| 24mm • 3/8 Circle | NFS-1 | FS-1 | C-14 | DS24 | C7 |
| 26mm • 3/8 Circle | NFS | FS | C-15 | DS26 | C8 |
| 30mm • 3/8 Circle | NFSL | FSL | C-16 | DS30 | C10 |
| 40mm • 3/8 Circle | NFSLX | FSLX | C-17 | DS40 | • |
| 13mm • 1/2 Circle | NM-1 | M-1 | C-21 | • | C26 |
| 15mm • 1/2 Circle | NX-8 | X-8 | • | • | • |
| 17mm • 3/8 Circle | NJ-1 | J-1 | C-26 | • | C20 |
| 23mm • 3/8 Circle | NX-1 | X-1 | C-23 | HS24 | • |
| 26mm • 1/2 Circle | NCP-2 | CP-2 | GS-10 | • | C21 |
| 37mm • 1/2 Circle | NCP-1 | CP-1 | GS-11 | • | C12 |



[Needles are to scale]

STRAIGHT CUTTING NEEDLES ▽

| NEEDLE | CP MEDICAL | ETHICON | COVIDIEN | SHARPOINT | LOOK |
|--------|------------|---------|----------|-----------|------|
| 60mm | NKS | KS | SC-2 | CS-1 | • |



TAPER POINT NEEDLES ○

| NEEDLE | CP MEDICAL | ETHICON | COVIDIEN | SHARPOINT | LOOK |
|-------------------|------------|---------|----------|-----------|------|
| 17mm • 3/8 Circle | NBB | BB | CV-15 | • | • |
| 17mm • 1/2 Circle | NRB-1 | RB-1 | CV-23 | HR17 | T28 |
| 22mm • 1/2 Circle | NCT-3 | CT-3 | GS-23 | • | • |
| 22mm • 1/2 Circle | NSH-1 | SH-1 | CV-25 | HR22 | T4 |
| 26mm • 1/2 Circle | NSH | SH | V-20 | HR26 | T |
| 26mm • 1/2 Circle | NMO-6 | MO-6 | HGS-22 | • | • |
| 26mm • 1/2 Circle | NCT-2 | CT-2 | GS-22 | • | T5 |
| 37mm • 1/2 Circle | NCT-1 | CT-1 | GS-21 | • | T8 |
| 40mm • 1/2 Circle | NCT | CT | GS-24 | • | T10 |



SOFTCUT NEEDLES ◎

| NEEDLE | CP MEDICAL | ETHICON | COVIDIEN | SHARPOINT | LOOK |
|-------------------|------------|---------|----------|-----------|------|
| 17mm • 3/8 Circle | NV-4 | V-4 | KV-15 | • | • |
| 26mm • 1/2 Circle | NV-7 | V-7 | KV-7 | • | • |



CONVENTIONAL CUTTING STERNUM NEEDLES △

| NEEDLE | CP MEDICAL | ETHICON | COVIDIEN | SHARPOINT | LOOK |
|-------------------|------------|---------|----------|-----------|------|
| 48mm • 1/2 Circle | NCCS | CCS | SCC | • | • |



[Needles are to scale]